

For Research Use Only

ADAMTSL3 Polyclonal antibody, PBS Only

Catalog Number: 26653-1-PBS



Basic Information

Catalog Number: 26653-1-PBS	GenBank Accession Number: BC128389	Purification Method: Antigen affinity purification
Concentration: 1 mg/ml	GeneID (NCBI): 57188	
Source: Rabbit	UNIPROT ID: P82987	
Isotype: IgG	Full Name: ADAMTS-like 3	
Immunogen Catalog Number: AG24488	Calculated MW: 1691 aa, 189 kDa	

Applications

Tested Applications:
IHC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

ADAMTSL3 (ADAMTS-Like Protein 3) is a member of the ADAMTSL family of matricellular proteins, which are involved in extracellular matrix (ECM) remodeling and regulation of various biological processes, including cardiac function, cartilage development, and inflammation. (PMID: 38324186)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 26653-1-AP (ADAMTSL3 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26653-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 26653-1-AP (ADAMTSL3 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26653-1-PBS in a different storage buffer formulation.