

For Research Use Only

SSPN Polyclonal antibody

Catalog Number: 26685-1-AP



Basic Information

Catalog Number:

26685-1-AP

Size:

1200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24511

GenBank Accession Number:

BC062299

GeneID (NCBI):

8082

UNIPROT ID:

Q14714

Full Name:

sarcospan (Kras oncogene-associated gene)

Calculated MW:

243 aa, 27 kDa

Observed MW:

15-20 kDa, 27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26685-1-AP (SSPN Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.