

For Research Use Only

FGF-7/KGF Polyclonal antibody, PBS Only

Catalog Number: 26698-1-PBS



Basic Information

Catalog Number:

26698-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24903

GenBank Accession Number:

BC010956

GeneID (NCBI):

2252

UNIPROT ID:

P21781

Full Name:

fibroblast growth factor 7
(keratinocyte growth factor)

Calculated MW:

22-28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 26698-1-AP (FGF7 Antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 26698-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 26698-1-AP (FGF7 Antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 26698-1-PBS in a different storage buffer formulation.