

For Research Use Only

FLVCR2 Polyclonal antibody

Catalog Number: 26704-1-AP **1 Publications**



Basic Information

Catalog Number: 26704-1-AP	GenBank Accession Number: BC019087	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 55640	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9UPI3	
Isotype: IgG	Full Name: feline leukemia virus subgroup C cellular receptor family, member 2	
Immunogen Catalog Number: AG24928	Calculated MW: 526 aa, 57 kDa	
	Observed MW: 65-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, rat liver
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiuhua Su	36313565	Front Cell Dev Biol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26704-1-AP (FLVCR2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.