

For Research Use Only

# LAMB3 Polyclonal antibody

Catalog Number: 26795-1-AP

Featured Product

14 Publications



## Basic Information

Catalog Number:

26795-1-AP

Concentration:

600 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25192

GenBank Accession Number:

BC075838

GeneID (NCBI):

3914

UNIPROT ID:

Q13751

Full Name:

laminin, beta 3

Calculated MW:

1172 aa, 130 kDa

Observed MW:

130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : A431 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xiaoyu Sui	36305754	Cancer Med	WB
Caiqun Huang	33179081	Mol. Med Rep	WB
Emma Plana	32605321	Int J Mol Sci	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 26795-1-AP (LAMB3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.