

For Research Use Only

MRPS7 Polyclonal antibody

Catalog Number: 26828-1-AP

3 Publications



Basic Information

Catalog Number:

26828-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25325

GenBank Accession Number:

BC000241

GeneID (NCBI):

51081

UNIPROT ID:

Q9Y2R9

Full Name:

mitochondrial ribosomal protein S7

Calculated MW:

28 kDa

Observed MW:

27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HepG2 cells, PC-3 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ying Wang	34707288	Nature	WB
Shabnam Bakhshalizadeh	37148394	Hum Genet	WB
Sumudu S C Amarasekera	37133451	Hum Mol Genet	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26828-1-AP (MRPS7 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.