

For Research Use Only

# PLB1 Polyclonal antibody, PBS Only

Catalog Number: 26835-1-PBS



## Basic Information

Catalog Number:

26835-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25240

GenBank Accession Number:

BC150633

GeneID (NCBI):

151056

UNIPROT ID:

Q6P1J6

Full Name:

phospholipase B1

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



DU 145 cells were subjected to SDS PAGE followed by western blot with 26835-1-AP (PLB1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate hyperplasia tissue slide using 26835-1-AP (PLB1 Antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate hyperplasia tissue slide using 26835-1-AP (PLB1 Antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.