

For Research Use Only

Bora Polyclonal antibody

Catalog Number: 26866-1-AP **2 Publications**



Basic Information

Catalog Number: 26866-1-AP	GenBank Accession Number: BC056876	Purification Method: Antigen affinity purification
Size: 850 µg/ml	GeneID (NCBI): 79866	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q6PGQ7	
Isotype: IgG	Full Name: chromosome 13 open reading frame 34	
Immunogen Catalog Number: AG19123	Calculated MW: 559 aa, 61 kDa	
	Observed MW: 80-85 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells, Jurkat cells
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Bin Yu	34313310	J Cell Sci	WB
Bin Yu	34060616	J Cell Sci	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 26866-1-AP (Bora antibody at dilution of 1:1000 incubated at room temperature for 1.5 hours.