

For Research Use Only

Bora Polyclonal antibody

Catalog Number: 26866-1-AP

2 Publications



Basic Information

Catalog Number:

26866-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19123

GenBank Accession Number:

BC056876

GeneID (NCBI):

79866

UNIPROT ID:

Q6PGQ7

Full Name:

chromosome 13 open reading frame

34

Calculated MW:

559 aa, 61 kDa

Observed MW:

80-85 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : K-562 cells, Jurkat cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Bin Yu	34313310	J Cell Sci	WB
Bin Yu	34060616	J Cell Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 26866-1-AP (Bora antibody at dilution of 1:1000 incubated at room temperature for 1.5 hours.