

For Research Use Only

YY2 Polyclonal antibody

Catalog Number: 26878-1-AP



Basic Information

Catalog Number: 26878-1-AP	GenBank Accession Number: BC137217	Purification Method: Antigen affinity purification
Size: 370 µg/ml	GeneID (NCBI): 404281	
Source: Rabbit	UNIPROT ID: O15391	
Isotype: IgG	Full Name: YY2 transcription factor	
Immunogen Catalog Number: AG25373	Calculated MW: 372 aa, 46 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data