

For Research Use Only

CNTROB Polyclonal antibody, PBS Only

Catalog Number: 26880-1-PBS



Basic Information

Catalog Number:

26880-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25454

GenBank Accession Number:

BC021134

GeneID (NCBI):

116840

UNIPROT ID:

Q8N137

Full Name:

centrobin, centrosomal BRCA2
interacting protein

Calculated MW:

903 aa, 101 kDa

Observed MW:

60-70 kDa, 100-130 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 26880-1-AP (CNTROB Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26880-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed hTERT-RPE1 cells using CNTROB antibody (26880-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 26880-1-PBS in a different storage buffer formulation.