

For Research Use Only

Nurim Polyclonal antibody, PBS Only

Catalog Number: 27009-1-PBS



Basic Information

Catalog Number:

27009-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25501

GenBank Accession Number:

BC039865

GeneID (NCBI):

11270

UNIPROT ID:

Q8IXM6

Full Name:

nurim (nuclear envelope membrane protein)

Calculated MW:

29 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

Nurim (nuclear envelope membrane protein) is also known as NRM and NRM29. Its primary cellular localization is the inner nuclear membrane (INM), where it is tightly associated with the nucleus. Nurim has emerged as a crucial regulator of apoptosis suppression.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27009-1-AP (Nurim antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27009-1-PBS in a different storage buffer formulation.