

For Research Use Only

# RNF113A Polyclonal antibody, PBS Only

Catalog Number: 27018-1-PBS

Featured Product



## Basic Information

Catalog Number:

27018-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25313

GenBank Accession Number:

BC000832

GeneID (NCBI):

7737

UNIPROT ID:

O15541

Full Name:

ring finger protein 113A

Calculated MW:

39 kDa

Observed MW:

40-50 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27018-1-AP (RNF113A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27018-1-PBS in a different storage buffer formulation.



WB result of RNF113A antibody (27018-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-RNF113A transfected MCF-7 cells. This data was developed using the same antibody clone with 27018-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 27018-1-AP (RNF113A antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 27018-1-PBS in a different storage buffer formulation.