

For Research Use Only

TAF5 Polyclonal antibody, PBS Only

Catalog Number: 27041-1-PBS



Basic Information

Catalog Number:

27041-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25742

GenBank Accession Number:

BC136348

GeneID (NCBI):

6877

UNIPROT ID:

Q15542

Full Name:

TAF5 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 100kDa

Calculated MW:

87 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27041-1-AP (TAF5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27041-1-PBS in a different storage buffer formulation.



WB result of TAF5 antibody (27041-1-AP; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-TAF5 transfected HEK-293 cells. This data was developed using the same antibody clone with 27041-1-PBS in a different storage buffer formulation.