

For Research Use Only

ZNF367 Polyclonal antibody

Catalog Number: 27046-1-AP



Basic Information

Catalog Number:

27046-1-AP

Size:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23263

GenBank Accession Number:

BC126141

GeneID (NCBI):

195828

UNIPROT ID:

Q7RTV3

Full Name:

zinc finger protein 367

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HL-60 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

