

For Research Use Only

SNAPC3 Polyclonal antibody

Catalog Number: 27062-1-AP



Basic Information

Catalog Number:

27062-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19851

GenBank Accession Number:

BC108743

GeneID (NCBI):

6619

UNIPROT ID:

Q92966

Full Name:

small nuclear RNA activating complex, polypeptide 3, 50kDa

Calculated MW:

411 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, HEK-293T cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27062-1-AP (SNAPC3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.