

For Research Use Only

SHCBP1L Polyclonal antibody

Catalog Number: 27108-1-AP



Basic Information

Catalog Number: 27108-1-AP	GenBank Accession Number: BC050305	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 81626	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9BZQ2	
Isotype: IgG	Full Name: chromosome 1 open reading frame 14	
Immunogen Catalog Number: AG25903	Calculated MW: 725 aa, 81 kDa	
	Observed MW: 81 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Species Specificity: human, mouse, rat	

Background Information

SHCBP1L, also known as C1orf14, is predominantly expressed in male germ cells. In association with HSPA2, it participates in the maintenance of spindle integrity during meiosis in male germ cells.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 27108-1-AP (SHCBP1L antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.