

For Research Use Only

PRMT5 Polyclonal antibody

Catalog Number: 27109-1-AP



Basic Information

Catalog Number:

27109-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25896

GenBank Accession Number:

BC025979

GeneID (NCBI):

10419

UNIPROT ID:

O14744

Full Name:

protein arginine methyltransferase 5

Calculated MW:

637 aa, 73 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HEK-293 cells, PC-12 cells, NIH/3T3 cells, Raji cells, HeLa cells, HepG2 cells, HL-60 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cells were subjected to SDS PAGE followed by western blot with 27109-1-AP (PRMT5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.