

For Research Use Only

SIVA1 Polyclonal antibody

Catalog Number: 27119-1-AP



Basic Information

Catalog Number:

27119-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25748

GenBank Accession Number:

BC034562

GeneID (NCBI):

10572

UNIPROT ID:

O15304

Full Name:

SIVA1, apoptosis-inducing factor

Calculated MW:

175 aa, 19 kDa

Observed MW:

19 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MOLT-4 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27119-1-AP (SIVA1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.