

For Research Use Only

# CCDC90B Polyclonal antibody

Catalog Number: 27126-1-AP



## Basic Information

Catalog Number:

27126-1-AP

Size:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25934

GenBank Accession Number:

BC014573

GeneID (NCBI):

60492

UNIPROT ID:

Q9GZT6

Full Name:

coiled-coil domain containing 90B

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : MCF-7 cells, human placenta, mouse testis tissue, rat spleen tissue

## Background Information

CCDC90B (Coiled-coil domain-containing protein 90B) is predicted to have a single-pass transmembrane domain located at its C-terminal. As two isoforms belong to the same protein family, CCDC90A/MCUR1 and CCDC90B show similar expression levels among tissues with few exceptions.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27126-1-AP (CCDC90B antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.