

For Research Use Only

PSRC1 Polyclonal antibody

Catalog Number: 27130-1-AP



Basic Information

Catalog Number:

27130-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26039

GenBank Accession Number:

BC056909

GeneID (NCBI):

84722

UNIPROT ID:

Q6PGN9

Full Name:

proline/serine-rich coiled-coil 1

Calculated MW:

39 kDa

Observed MW:

39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells, HEK-293T cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 27130-1-AP (PSRC1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.