

For Research Use Only

ABCB8 Polyclonal antibody

Catalog Number: 27158-1-AP



Basic Information

Catalog Number:

27158-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25935

GenBank Accession Number:

BC141814

GeneID (NCBI):

11194

UNIPROT ID:

Q9NUT2

Full Name:

ATP-binding cassette, sub-family B (MDR/TAP), member 8

Calculated MW:

735 aa, 80 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

