

For Research Use Only

ZDHHC8 Polyclonal antibody

Catalog Number: 27179-1-AP **2 Publications**



Basic Information

Catalog Number: 27179-1-AP	GenBank Accession Number: BC053544	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 29801	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q9ULC8	
Isotype: IgG	Full Name: zinc finger, DHHC-type containing 8	
Immunogen Catalog Number: AG26028	Calculated MW: 83 kDa	
	Observed MW: 70-85 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, K-562 cells, rat brain tissue
Cited Applications: WB, IHC	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhangjie Wang	38211651	Cancer Lett	WB,IHC
Ryosuke Shimai	37544935	Sci Rep	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27179-1-AP (ZDHC8 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.