

For Research Use Only

CCDC137 Polyclonal antibody

Catalog Number: 27201-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

27201-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23942

GenBank Accession Number:

BC009369

GeneID (NCBI):

339230

UNIPROT ID:

Q6PK04

Full Name:

coiled-coil domain containing 137

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : PC-3 cells, HEK-293 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shuang Tao	37542342	J Exp Clin Cancer Res	WB,IHC,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 27201-1-AP (CCDC137 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.