

For Research Use Only

# LMBRD1 Polyclonal antibody, PBS Only

Catalog Number: 27203-1-PBS



## Basic Information

Catalog Number:

27203-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24478

GenBank Accession Number:

BC010360

GeneID (NCBI):

55788

UNIPROT ID:

Q9NUN5

Full Name:

LMBR1 domain containing 1

Calculated MW:

540 aa, 61 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 27203-1-AP (LMBRD1 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 27203-1-AP (LMBRD1 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27203-1-PBS in a different storage buffer formulation.