

For Research Use Only

CCDC114 Polyclonal antibody

Catalog Number: 27211-1-AP



Basic Information

Catalog Number:

27211-1-AP

Size:

1300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26053

GenBank Accession Number:

NM_144577

GeneID (NCBI):

93233

UNIPROT ID:

Q96M63

Full Name:

coiled-coil domain containing 114

Calculated MW:

75 kDa

Observed MW:

65kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27211-1-AP (CCDC114 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.