

For Research Use Only

CCDC114 Polyclonal antibody

Catalog Number: 27211-1-AP



Basic Information

Catalog Number: 27211-1-AP	GenBank Accession Number: NM_144577	Purification Method: Antigen affinity purification
Size: 1300 µg/ml	GeneID (NCBI): 93233	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q96M63	
Isotype: IgG	Full Name: coiled-coil domain containing 114	
Immunogen Catalog Number: AG26053	Calculated MW: 75 kDa	
	Observed MW: 65kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27211-1-AP (CCDC114 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.