

For Research Use Only

RAB5C Polyclonal antibody

Catalog Number: 27219-1-AP **1 Publications**



Basic Information

Catalog Number: 27219-1-AP	GenBank Accession Number: BC106039	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 5878	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: P51148	
Isotype: IgG	Full Name: RAB5C, member RAS oncogene family	
Immunogen Catalog Number: AG25849	Calculated MW: 216 aa, 23 kDa	
	Observed MW: 23 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, U-937 cells
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jun Yan	34367131	Front Immunol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 27219-1-AP (RAB5C antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.