

For Research Use Only

AS3MT Polyclonal antibody

Catalog Number: 27270-1-AP



Basic Information

Catalog Number:

27270-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26194

GenBank Accession Number:

NM_020682

GeneID (NCBI):

57412

Full Name:

arsenic (+3 oxidation state)

methyltransferase

Calculated MW:

42 kDa

Observed MW:

39-41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 27270-1-AP (AS3MT Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.