

For Research Use Only

C14orf153 Polyclonal antibody

Catalog Number: 27300-1-AP **1 Publications**



Basic Information

Catalog Number: 27300-1-AP	GenBank Accession Number: BC007412	Purification Method: Antigen affinity purification
Size: 270 µg/ml	GeneID (NCBI): 84334	
Source: Rabbit	UNIPROT ID: Q961L0	
Isotype: IgG	Full Name: chromosome 14 open reading frame 153	
Immunogen Catalog Number: AG23870		

Applications

Tested Applications:
IF/ICC, ELISA

Cited Applications:
WB, IF

Species Specificity:
human

Cited Species:
human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Alba Signes	30552096	EMBO Mol Med	IF, WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data