

For Research Use Only

ETV1 Polyclonal antibody

Catalog Number: 27381-1-AP **1 Publications**



Basic Information

Catalog Number:

27381-1-AP

Size:

470 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26499

GenBank Accession Number:

BC098403

GeneID (NCBI):

2115

UNIPROT ID:

P50549

Full Name:

ets variant 1

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB: MCF-7 cells, mouse testis tissue, mouse brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yangsong Zuo	32308483	Cancer Manag Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 27381-1-AP (ETV1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.