

For Research Use Only

# WDR91 Polyclonal antibody

Catalog Number: 27424-1-AP



## Basic Information

Catalog Number:

27424-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26697

GenBank Accession Number:

BC017246

GeneID (NCBI):

29062

UNIPROT ID:

A4D1P6

Full Name:

WD repeat domain 91

Observed MW:

42 kDa, 83 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, Jurkat cells, SH-SY5Y cells, U2OS cells

## Background Information

WDR91 (1-405 aa) can be detected 42 kDa (PMID:28860274).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27424-1-AP (WDR91 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.