

For Research Use Only

SH3BP5L Polyclonal antibody

Catalog Number: 27446-1-AP



Basic Information

Catalog Number: 27446-1-AP	GenBank Accession Number: BC017254	Purification Method: Antigen affinity purification
Concentration: 500 ug/ml	GeneID (NCBI): 80851	Recommended Dilutions: WB: 1:2000-1:16000
Source: Rabbit	UNIPROT ID: Q7L8J4	
Isotype: IgG	Full Name: SH3-binding domain protein 5-like	
Immunogen Catalog Number: AG23990	Observed MW: 65-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : 769-P cells, PC-3 cells, RKO cells, THP-1 cells
Species Specificity: human	

Background Information

SH3BP5L (SH3 Binding Domain Protein 5 Like) is a protein-coding gene that encodes a protein with a conserved N-terminus and functions as a guanine nucleotide exchange factor (GEF) for Rab11 GTPases.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27446-1-AP (SH3BP5L antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.