

For Research Use Only

# HES1 Polyclonal antibody

Catalog Number: 27542-1-AP

4 Publications



## Basic Information

Catalog Number:

27542-1-AP

Size:

370 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26393

GenBank Accession Number:

BC039152

GeneID (NCBI):

3280

UNIPROT ID:

Q14469

Full Name:

hairy and enhancer of split 1,  
(Drosophila)

Calculated MW:

30 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human, mouse, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Tianrong Zhang	31391046	Lipids Health Dis	WB
Xueying Jin	38703204	Cell Mol Life Sci	WB
Menghuan Wang	38134973	J Nutr Biochem	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

