

For Research Use Only

FAM53C Polyclonal antibody

Catalog Number: 27612-1-AP



Basic Information

Catalog Number:

27612-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14219

GenBank Accession Number:

BC000259

GeneID (NCBI):

51307

UNIPROT ID:

Q9NYF3

Full Name:

family with sequence similarity 53, member C

Calculated MW:

392 aa, 43 kDa

Observed MW:

43-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : K-562 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells lysates were subjected to SDS PAGE followed by western blot with 27612-1-AP (FAM53C antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.