

For Research Use Only

# WDR8 Polyclonal antibody

Catalog Number: 27628-1-AP



## Basic Information

Catalog Number:	27628-1-AP	GenBank Accession Number:	BC002611	Purification Method:	Antigen affinity purification
Size:	500 µg/ml	GeneID (NCBI):	49856	Recommended Dilutions:	IHC 1:50-1:500
Source:	Rabbit	UNIPROT ID:	Q9P2S5		
Isotype:	IgG	Full Name:	WD repeat domain 8		
Immunogen Catalog Number:	AG26267	Calculated MW:	460 aa, 52 kDa		

## Applications

Tested Applications:	IHC, ELISA	Positive Controls:	
Species Specificity:	Human	IHC :	human heart tissue, human brain tissue
<b>Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0</b>			

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 27628-1-AP (WDR8 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 27628-1-AP (WDR8 antibody) at dilution of 1:200 (under 40x lens).