

For Research Use Only

WDR8 Polyclonal antibody

Catalog Number: 27628-1-AP



Basic Information

Catalog Number: 27628-1-AP	GenBank Accession Number: BC002611	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 49856	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q9P255	
Isotype: IgG	Full Name: WD repeat domain 8	
Immunogen Catalog Number: AG26267	Calculated MW: 460 aa, 52 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
Human

Positive Controls:

IHC : human heart tissue, human brain tissue

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 27628-1-AP (WDR8 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 27628-1-AP (WDR8 antibody) at dilution of 1:200 (under 40x lens).