

For Research Use Only

COX14 Polyclonal antibody

Catalog Number: 27712-1-AP



Basic Information

Catalog Number: 27712-1-AP	GenBank Accession Number: BC007849	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 84987	Recommended Dilutions: WB: 1:500-1:1000 IF/ICC: 1:50-1:500
Isotype: IgG	UNIPROT ID: Q96I36	
Immunogen Catalog Number: AG24036	Full Name: Cytochrome c oxidase assembly protein COX14	
	Observed MW: 7-12 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HEK-293 cells, HepG2 cells IF/ICC : HepG2 cells,
Species Specificity: human	

Background Information

COX14 (Cytochrome c Oxidase Assembly Factor 14) is a protein that functions as a crucial assembly factor for mitochondrial complex IV (cytochrome c oxidase, COX), the terminal enzyme of the electron transport chain.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27712-1-AP (COX14 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using COX14 antibody (27712-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red).