

For Research Use Only

# ABCB4 Polyclonal antibody

Catalog Number: 27726-1-AP

1 Publications



## Basic Information

Catalog Number:

27726-1-AP

Concentration:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26840

GenBank Accession Number:

BC042531

GeneID (NCBI):

5244

UNIPROT ID:

P21439

Full Name:

ATP-binding cassette, sub-family B (MDR/TAP), member 4

Calculated MW:

1286 aa, 142 kDa

Observed MW:

130-140 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yuhan Wang	39578442	Cell Death Dis	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27726-1-AP (ABCB4 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.