

For Research Use Only

FAM91A1 Polyclonal antibody

Catalog Number: 27738-1-AP

1 Publications



Basic Information

Catalog Number:

27738-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26902

GenBank Accession Number:

BC030520

GeneID (NCBI):

157769

UNIPROT ID:

Q658Y4

Full Name:

family with sequence similarity 91, member A1

Observed MW:

94 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IP

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : A549 cells,

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------------|-----------|---------|-------------|
| Huaqing Deng | 39013469 | Cell | IP |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 and HeLa cells were subjected to SDS PAGE followed by western blot with 27738-1-AP (FAM91A1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.