

For Research Use Only

CCDC147 Polyclonal antibody, PBS Only

Catalog Number: 27739-1-PBS



Basic Information

Catalog Number: 27739-1-PBS	GenBank Accession Number: BC109126	Purification Method: Antigen affinity purification
Size: 1 mg/ml	GeneID (NCBI): 159686	
Source: Rabbit	UNIPROT ID: Q5T655	
Isotype: IgG	Full Name: coiled-coil domain containing 147	
Immunogen Catalog Number: AG26943	Calculated MW: 872 aa, 103 kDa	
	Observed MW: 103 kDa	

Applications

Tested Applications:
IP, Indirect ELISA

Species Specificity:
human, mouse

Background Information

CCDC147 has an essential role in the assembly and organization of the sperm flagellar axoneme (PMID: 32791035).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-CCDC147 (IP:27739-1-AP, 4ug; Detection:27739-1-AP 1:1000) with mouse testis tissue lysate 1600 ug. This data was developed using the same antibody clone with 27739-1-PBS in a different storage buffer formulation.