

For Research Use Only

DDIAS Polyclonal antibody, PBS Only

Catalog Number: 27754-1-PBS



Basic Information

Catalog Number:

27754-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26829

GenBank Accession Number:

BC039268

GeneID (NCBI):

220042

UNIPROT ID:

Q8IXT1

Full Name:

chromosome 11 open reading frame

82

Observed MW:

120 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

DNA damage-induced apoptosis suppressor (DDIAS) is an oncogenic protein that is highly expressed in several cancers, including colorectal, lung, breast, and hepatocellular carcinoma (HCC), and stimulates cancer cell proliferation and cell cycle progression. DDIAS protects cancer cells from DNA damage reagents and tumor necrosis factor-associated apoptosis-inducing ligands (TRAILs) and positively regulates cancer cell invasion by stabilizing the β -catenin. β -catenin positively regulates cancer cell invasion. Therefore, DDIAS is considered as a potential therapeutic target for malignant lung cancer.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27754-1-AP (C11orf82 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27754-1-PBS in a different storage buffer formulation.