

For Research Use Only

C17orf80 Polyclonal antibody

Catalog Number: 27762-1-AP



Basic Information

Catalog Number: 27762-1-AP	GenBank Accession Number: BC005005	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 55028	Recommended Dilutions: WB 1:1000-1:8000
Source: Rabbit	UNIPROT ID: Q9BSJ5	
Isotype: IgG	Full Name: chromosome 17 open reading frame 80	
Immunogen Catalog Number: AG26903	Calculated MW: 609 aa, 66 kDa	
	Observed MW: 66 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 27762-1-AP (C17orf80 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.