

For Research Use Only

C17orf80 Polyclonal antibody

Catalog Number: 27762-1-AP

1 Publications



Basic Information

Catalog Number:

27762-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26903

GenBank Accession Number:

BC005005

GeneID (NCBI):

55028

UNIPROT ID:

Q9BSJ5

Full Name:

chromosome 17 open reading frame
80

Calculated MW:

609 aa, 66 kDa

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : HEK-293 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hao Wu	37676315	J Cell Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 27762-1-AP (C17orf80 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.