

For Research Use Only

# MTNAP1 Polyclonal antibody, PBS Only

Catalog Number: 27762-1-PBS



## Basic Information

Catalog Number:

27762-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26903

GenBank Accession Number:

BC005005

GeneID (NCBI):

55028

UNIPROT ID:

Q9BSJ5

Full Name:

chromosome 17 open reading frame

80

Calculated MW:

609 aa, 66 kDa

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 27762-1-AP (C17orf80 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27762-1-PBS in a different storage buffer formulation.