

For Research Use Only

RRAD Polyclonal antibody

Catalog Number: 27763-1-AP



Basic Information

Catalog Number:

27763-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26869

GenBank Accession Number:

BC057815

GeneID (NCBI):

6236

UNIPROT ID:

P55042

Full Name:

Ras-related associated with diabetes

Calculated MW:

33 kDa

Observed MW:

29-34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : mouse heart tissue, rat heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27763-1-AP (RRAD antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.