

For Research Use Only

RIC3 Polyclonal antibody

Catalog Number: 27776-1-AP

2 Publications



Basic Information

Catalog Number:

27776-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26740

GenBank Accession Number:

BC022455

GeneID (NCBI):

79608

UNIPROT ID:

Q7Z5B4

Full Name:

resistance to inhibitors of
cholinesterase 3 homolog (C.
elegans)

Calculated MW:

368 aa, 41 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

Human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse heart tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Cai-Mei Zheng	33114531	Toxics	IHC
Dandan Wang	39639217	BMC Cancer	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 27776-1-AP (RIC3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.