

For Research Use Only

C9orf64 Polyclonal antibody, PBS Only

Catalog Number: 27797-1-PBS

Featured Product



Basic Information

Catalog Number:

27797-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23069

GenBank Accession Number:

BC004407

GeneID (NCBI):

84267

UNIPROT ID:

Q5T6V5

Full Name:

chromosome 9 open reading frame 64

Calculated MW:

341 aa, 39 kDa

Observed MW:

34-39 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of C9orf64 antibody (27797-1-AP; 1:20000; incubated at room temperature for 1.5 hours) with sh-Control and sh-C9orf64 transfected HeLa cells. This data was developed using the same antibody clone with 27797-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human intrahepatic cholangiocarcinoma tissue slide using 27797-1-AP (C9orf64 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 27797-1-PBS in a different storage buffer formulation.