

For Research Use Only

ABP1 Polyclonal antibody

Catalog Number: 27814-1-AP



Basic Information

Catalog Number:

27814-1-AP

Size:

750 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26057

GenBank Accession Number:

BC014093

GeneID (NCBI):

26

UNIPROT ID:

P19801

Full Name:

amiloride binding protein 1 (amine oxidase (copper-containing))

Calculated MW:

751 aa, 85 kDa

Observed MW:

85-87 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human placenta tissue, human colon cancer tissue, human kidney tissue, human stomach cancer tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 27814-1-AP (ABP1 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 27814-1-AP (ABP1 antibody) at dilution of 1:200 (under 40x lens).