

For Research Use Only

VSIG4 Polyclonal antibody

Catalog Number:27829-1-AP



Basic Information

Catalog Number:

27829-1-AP

Concentration:

700 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27214

GenBank Accession Number:

BC010525

GeneID (NCBI):

11326

UNIPROT ID:

Q9Y279

Full Name:

V-set and immunoglobulin domain containing 4

Calculated MW:

399 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : TT cells, U-87 MG cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27829-1-AP (V5IG4 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.